

# stayClean

Practical, single-use cover drape.



Order No. 30340014/Order No. 30160001

## Take efficiency a step further

To get the operating table up and running again as quickly as possible, we recommend to use the “stayClean”-surgical drapes by UFSK-OSYS.

Available in two different-sized custom fits, they provide a perfectly aseptic liner for every patient. They are impervious to fluids and minimized undesirable soiling of the upholstery. They protect the upholstery material of your treatment chair actively against aggressive fluids and have the treatment chair up and running again much faster. In addition, these surgical drapes keeps multi-resistant germs safely away from your patients.




Order No. 30340014




Order No. 30160001

## EQUIPMENT



1. Place the chair cover with the marked arrow facing upwards in the middle of the chair.



2. Unfold first to the left, then repeat from inside to outside.



3. Place the hood over the headrest.



4. Secure the straps onto the corresponding attachment points on both sides. Then move the cover downwards over the foot section.

## ADVANTAGES



### Advantages designed to increase your efficiency

- Easy replacement of used drapes
- Secure fastening thanks to cap, straps, and pre-punched holes
- Enormous time savings
- Low purchasing costs
- The new dispenser boxes are space-saving and fit on every surgery warehouse shelf



### Perfect protection of the equipment

- Protection against corrosion
- 100% waterproof rear side, absorptive front side
- Minimized soiling of the upholstery
- No sliding thanks to the head pocket



### Perfect protection of the patients

- Creates a barrier against multi-resistant germs
- Protects against cross-contamination

## SUPPLY

### Available in two sizes

- Surgical drape, short  
(covering the head and backrest elements)  
100 x 73 cm, 240 units (4 x 60 units), **Order No. 30160001**
- Surgical drape, long  
(covering the headrest through footrest element)  
200 x 73 cm, 160 units (4 x 40 units), **Order No. 30340014**

### Simply select the ordering mode and supply rhythm of your choice

- Complete delivery or partial delivery
- Partial deliveries are optionally available either as 'on-time deliveries' (deliveries made in accordance with delivery dates agreed in advance) or as 'deliveries on call' by phone, fax, e-mail, or via [www.ufsk-osys.com](http://www.ufsk-osys.com)
- Minimum unit available on-call: 1 carton containing 4 dispensers

## 4 in 1 box - for more organization in your OR

